
8/95-Sf-

*Literatur:*

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Fig. 1

Bibliography:

Physics Experiments, Volume 3 (599 942)
 New Physics Leaflets for Colleges and Universities,
 Volume 1 (599 952)
 Manual 'X-Ray Apparatus 42 kV' (554 911, in German)
 Quantum physics with the X-ray apparatus 42 kV (554 982)

1 Safety notes

Protect the fragile mica end window from mechanical damage, as this would make the end-window counter unusable.

Remove the protective cap only during measurement. Remove and replace the cap with care; do not twist the cap!

Do not touch the mica end-window.

Only store the counter with protective cap on.

Max. permissible operating voltage: 550 V d.c.
 (risk of damaging the counter by self-maintained discharge at voltages higher than 550 V).

2 Description, technical data

Self-quenching Geiger-Müller counter in plastic housing, with mica end-window;

sensitive to γ -rays, X-rays*) and high-energy**) β -rays;

mica end-window covered with protective cap;

with coaxial plug for counter-tube cable (559 06 or 07) to connect pulse counters with built-in voltage source for the counter.

*) The X-ray apparatuses 554 61 and 554 90 are already prepared for the installation of the end-window counter for β -, γ - and X-rays (559 05) as required for the experiments.

**) Unsuitable for the soft β -rays from cylindrical samples (from 559 35) artificially activated in the shielded source of neutrons (559 30, no longer available); counter-tube housing does not allow soft β -rays to enter the counter.

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
Gas filling: halogen
 Mean operating voltage: 450 V
 Plateau width: approx. 200 V
 Relative plateau slope: 0.05 %/V
 Dead time: approx. 100 µs
 Service life: 10¹⁰ pulses

Background in the plateau: approx. 0.2 pulses/s (with shielding of 50 mm Pb and 3 mm Al)
 approx. 1 %

Responsivity to γ-rays: 9 mm
 Diameter of mica window: 7.5 cm x 2.5 cm
 Dimensions:

3 Use

Observe Safety Notes (Section 1)!

Additionally required:
 Counter-tube cable, 25 cm or 100 cm 559 06 or 07

Pulse counter and voltage source (approx. 450 V DC); see Fig. 5
 Place the counter into the beam of radiation as required for the experiment (see examples in Figs. 2 to 4).

Connect the counter via counter-tube cable (559 06 or 07) to a pulse counter with built-in voltage source (approx. 450 V DC), see Fig. 5.

Fig. 2

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Range of β radiation; absorption measurements for β and γ radiation; holder for preparation in holder with absorber foils (559 18)


Fig. 3

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Deflection of β-rays by a magnetic field; rotatable fixture of counter in swivelling clamp (559 23)

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In experiments with low count rates (weak radioactive preparations, backscattering, absorption measurements) take the background radiation into account.

Indication of an overloaded counter not capable of processing the excessive number of particles entering the counter per unit of time: reduction of the measured pulse rate at increasing or constant radiation. Remedy: Reduce penetrating radiation (e. g. by increasing the distance by using a shield or absorber foils from 559 18).

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



Recommended pulse counters and readouts or counting devices, with built-in voltage source for powering the counter